

Included With Your Purchase

- 1 x Blaux Thermometer
- 1 x Instruction Manual

Normal Temperature Range By Age

Age	Celsius	Fahrenheit
0 - 2	36.4°- 38.0°	97.5°- 100.4°
3 - 10	36.1 °- 37.8°	97.0° - 100.0°
11 - 65	35.9°- 37.6°	96.6°- 99.7°
Over 65	35.8°- 37.5°	96.4°- 99.5°

Troubleshooting

Issue	Cause	Solution
Device Not Powering On	Dead battery	Replace the battery
	Batteries installed incorrectly	Ensure the batteries are installed correctly by making sure the Positive (+) and Negative(-) terminals are lined up.
Battery Symbol Displayed	Low battery	Replace the batteries
Incorrect Temperature Displayed	Device too far away from object when measurement taken	Ensure you are at the correct distance when taking your measurement 3-5cm (1-2")
	Device too cold	If you have left your thermometer in a cold place, it should be left at the current room temperature for 30 minutes
	If the result is far too high, the temperature head could be damaged or could have been exposed to direct sunlight	Contact manufacturer if damaged. Never leave the device under direct sunlight.

BLAUX


About Your Blaux Thermometer

The Blaux Thermometer is a non-contact infrared thermometer that easily and instantly measures internal body temperature using specialized infrared receiver technology.

To use it, simply point the sensor at your subject's forehead, roughly 3-5cm (1-2") away. Then, press the measurement key "SCAN" and your thermometer will instantly and accurately display the temperature on the screen.

More information on how to use all of these great features is listed in the "Basic Functions" section.

Structure and Function Key



The Science Behind Blaux Thermometer

All objects that have a temperature greater than 0°C (32°F) give off a small amount of infrared radiation.

For instance, a healthy temperature range of 36°C or 37°C will give off an infrared wavelength of 9-13µm.

Blaux Thermometer's infrared thermometer is able to measure and detect those wavelengths to give you accurate readings of temperature without the need to touch your patient.

Power Supply

Blaux Thermometer only needs 2 AAA batteries to operate.

*Battery not included

Working Specifications

- Temperature Units: Displays both°C and °F
- Measuring Distance: 3-5cm (1-2")
- Memory: 32 slots of temperature memory
- Detection Range: body temperature - 35.5°C-42.9°C; object - 0-100°C
- Accuracy± 0.2°C / ±0.4°F

Product Features

Safety

- Infrared technology does not produce harmful radiation
- 100% no-contact readings prevent the spread of infection

Ease of Use

- Ergonomic design fits easily in your hand
- One-click measurement
- Non-invasive (works well with sleeping children)

Quick and Accurate Measurement

- Measurements take only a fraction of a second
- Advanced infrared technology measures microscopic wavelengths for best accuracy

Special Design

- IR sensor lever can be rotated 90 degrees
- There is a tripod screw adapter at the bottom
- Slim and pocket-table design

Basic Functions

Powering On/Off	Press the "SCAN" key to power the device on.The device will automatically darken the display and shut down after 20 seconds of not being used.
Switching Between °C and °F	To switch between Celsius and Fahrenheit, press and hold the "C/F" button.
Taking a Measurement	When the device is on, point it at your subject and press the "SCAN" key, the temperature will be shown on the LCD screen. LCD Screen Lights: Green = Safe, Orange = Mild Fever, Red = Fever.
High Temperature Display (for human)	If the measured temperature exceeds 42°C (107.6°F), the device will simply display "HI" as the temperature.
Low Temperature Display (for human)	If the measured temperature is below 32°C (90°F), the device will simply display "LO" as the temperature.
Testing Scope	Body Temp: 35.5°C-42.9°C / 95.9°F-109°F Object Surface Temp: 0°C-100°C / 32°F-212°F
Memory Storage	To cycle through the previous readings, press the "▲" button. A number will be displayed on the screen representing previous measurements. For instance, "1" will represent the most recent recorded temperature, and "2" will represent the second most recent. After the memory number is displayed for one second, the temperature will be displayed, and you can cycle through the memory, which stores 32 readings. To clear the memory, press and hold the "▲" button until "CLR" appears. All memory slots will now be erased.
Low Battery	If your battery is running low when you try to start your device, "... " will be displayed on the screen. The battery symbol on the lower right corner will flash for two seconds, and the device will power off. At this point, you will need to swap in 2 new AAA batteries.

Have a question or concern?

Our customer support team is happy to help!

☎ US: 855 756 5448 / DE: 0800 6663955 / UK: 033081 80842 /
FR: 0185 16 9233 / AUS: (02) 5133 5672

✉ support@blaux.com 🌐 www.blaux.com