




IMPORTANT SAFETY INFORMATION

To prevent personal injury or damage to your KoreSurge™, read the following guidelines carefully to ensure proper usage:

1. If you are susceptible to heart rhythm disturbances, do not use KoreSurge™ across your chest.
2. If you are currently using other medical or physical treatment for muscle pain, consult a physician before using KoreSurge.
3. If while using KoreSurge™ your pain becomes more severe, stop using the device immediately.
4. Do not use KoreSurge™ in the bath or shower, or during any activity in which electrical stimulation can put you at risk of injury.

Warning: This device is not a toy. Keep away from children.




TIPS FOR USE:

- Start by placing KoreSurge™ on a sensitive or knotted spot and gradually increasing the vibration level to your comfort.
- Keep KoreSurge™ centered on your muscle tissue and avoid rolling over tendons, joints, or bony structures.
- If you find an area that's extra sensitive, focus KoreSurge™ on the muscles surround the spot rather than attacking it directly.
- Don't linger on one spot for too long. Spend 60 seconds or less on each sore area and then move on.

CHARGING THE BATTERY:

1. Do not use KoreSurge™ or activate vibration during charging.
2. Plug the charging cable into a wall outlet and inset the micro-usb plug into the charging port on the control panel.
3. If the power switch is in the ON position, the LED indicators will light up to show the charge level.
4. When fully charged, KoreSurge™ has a runtime of 20 hours.



NOTE: adaptor is not included in this package.


HOW TO USE:

1. Slide the power switch into the ON position.
2. Use the center power button to activate vibration and cycle through the different levels.
 - a. Press once to activate Level 1 – 900rpm
 - b. Press twice to activate Level 2 – 1800rpm
 - c. Press three times to activate Level 3 – 2800rpm
 - d. Press four times to activate Level 4 – 3800rpm
 - e. Press five times to activate Recovery Mode - auto switch from Mode 1 to Mode 4.
 - f. Press six times to return to Mode 1.

- Place KoreSurge™ on the floor and lower the sore or tight area of your body on top of it.
- Roll back and forth with KoreSurge™, slowly applying more pressure as needed.
- KoreSurge™ will automatically shut down after 15 minutes.

WORKOUT DEMONSTRATIONS:

Hamstrings




Outer Thighs




Knees




Inner Thighs



Back




Quads



SPECIFICATIONS:

Package Includes	KoreSurge™ foam roller, charging cable, and Instruction manual.
Product Size	Ø 14 x 33cm
Product Weight	1.35kg
Hardness	50 +/- 5 A
Material	Polypropylene, EVA
Massage Modes	Level 1 - 900rpm, Level 2 - 1800rpm, Level 3 - 2800rpm, Level 4 - 3800rpm, Level 5 - Recovery Mode
Battery Runtime	20 hours
Battery Capacity	4000mAh
Maximum Loading	120 kg
Input Port	MicroUSB
Input Voltage	5V1A

CE The device complies with the European Directives 2014/30/EU.

 At the end of its service life, this product may not be disposed of in the normal household waste, but rather must be disposed of at a collection point for recycling electrical and electronic devices. The materials are recyclable according to their labelling. You make an important contribution to protecting our environment by reusing, recycling, or utilizing old devices in other ways. Please ask the municipal administration where the appropriate disposal facility is located.

Subject to change.